



Centrella® Smart+ Bed:


**CONTACT-FREE,
CONTINUOUS
MONITORING**



Hillrom™

Contact-free, Continuous Monitoring

- ✓ Identify signs of patient deterioration
- ✓ Initiate an early intervention
- ✓ Promote patient safety



Promoting a Culture of Patient Safety



These organizations and others like them have recognized a simple fact:

TOO MANY PATIENTS ARE DYING FROM PREVENTABLE CAUSES.

2019 Top 10 Patient Safety Concerns¹

ECRI Institute

1. Diagnostic Stewardship and Test Result Management Using EHRs
2. Antimicrobial Stewardship in Physician Practices and Aging Services
3. Burnout and Its Impact on Patient Safety
4. Patient Safety Concerns Involving Mobile Health
5. Reducing Discomfort with Behavioral Health
6. Detecting Changes in a Patient's Condition
7. Developing and Maintaining Skills
8. Early Recognition of Sepsis Across the Continuum
9. Infections from Peripherally Inserted IV Lines
10. Standardizing Safety Efforts Across Large Health Systems

2019 Top 10 Patient Safety Concerns¹



1. Diagnostic Stewardship and Test Result Management Using EHRs
2. Antimicrobial Stewardship in Physician Practices and Aging Services
3. Burnout and Its Impact on Patient Safety
4. Patient Safety Concerns Involving Mobile Health
5. Reducing Discomfort with Behavioral Health
- 6. DETECTING CHANGES IN A PATIENT'S CONDITION**
7. Developing and Maintaining Skills
- 8. EARLY RECOGNITION OF SEPSIS ACROSS THE CONTINUUM**
9. Infections from Peripherally Inserted IV Lines
10. Standardizing Safety Efforts Across Large Health Systems

Critical Events Drive Patient Safety Outcomes



**UNRECOGNIZED PATIENT DETERIORATION
THAT CAN LEAD TO MORTALITY**

17%

MAY AFFECT AS MANY AS 17% OF HOSPITAL ADMISSIONS.²

Patient Deterioration

A patient moves from one clinical state to a worse clinical state.¹²

Increasing their individual risk of

MORBIDITY → PROTRACTED HOSPITAL STAY → DISABILITY → DEATH

FAILURE TO RESCUE¹³

Heart failure

Electrolyte abnormalities

Sepsis

Ischemia

DVT/PE


Respiratory insufficiency

Sepsis



#1 cause of death in U.S. hospitals⁸

35% OF ALL DEATHS IN HOSPITALS⁸



- Risk can be reduced by quickly identifying and managing infections.⁸
- Mortality increases **8%** for every hour that treatment is delayed.⁸
- **\$38k**, Median hospital cost to treat Hospital Acquired Severe Sepsis.¹⁰

Opioid Induced Respiratory Depression



57%

of medical patients were prescribed opioids, sedatives, or both.³

1/3

of Code Blue arrests are from respiratory depression.⁴


\$27k

increased treatment cost for opioid patient that suffers arrest.⁵

7.57 DAYS

increased length of stay of opioid patient suffering arrest.⁵

Patient Deterioration



6-8

POTENTIAL HOURS OF WARNING SIGNS PRIOR TO EVENT.²

The Med-Surg Environment is Challenging



INCREASING PATIENTS, ACUITY AND COMPLEXITY

STAGNANT TECHNOLOGY & INNOVATIONS

COST CONSTRAINTS


HIGH PATIENT-TO-NURSE RATIOS

STAFF TURNOVER, RETENTION AND TRAINING




Leading Indicators of Deterioration

Respiratory rate and heart rate are the most important predictor variables of deterioration.⁶



Two Paths of Patient Deterioration

CONTRIBUTING CONDITIONS



PATIENT
DETERIORATION

RECOGNIZED

RN Intervention • MD Intervention • Rapid Response

UNRECOGNIZED

Code Blue • ICU Transfer/Readmit • Mortality

A Real-Life Impact:⁹ Opioid Induced Respiratory Failure



John LaChance

- Rotator Cuff surgical patient
- Diagnosed with sleep apnea

Thursday, March 15, 2007

John underwent his second routine rotator cuff repair surgery.

10 AM

12 PM

2 PM

During the afternoon, John began to deteriorate. Symptoms included fever and extreme vomiting.

Within a half hour, John was comatose and never spoke again.

5 PM

5:30 PM

Immediately following surgery, John was doing well. His pain was managed with a shoulder block and Morphine through a PCA pump.

By the early evening, his shoulder block had worn off. John was taken off Morphine and prescribed a high dose of Dilaudid.




5:30 AM

Friday,
March 16, 2007

By 5:30am, John was dead.

Continuous Monitoring in an Inpatient Medical-Surgical Unit: A Controlled Clinical Trial

Harvey Brown, MD²



“Results may support the hypothesis that continuous monitoring leads to earlier recognition of patient deterioration.”

Identifying Patient Deterioration Early Using Contact-Free, Continuous Monitoring on an Inpatient Medical – Surgical Unit

Annot Health, accepted to be presented at the Institute for Healthcare Improvement Conference, December 2019

POSTER HIGHLIGHTS

- 26-bed inpatient medical-surgical unit
- 7 cases in which contact-free continuous monitoring was used to **identify patient deterioration** and **intervene earlier**

CASE TYPES

- Airway management
- Alcohol withdrawal management
- Opioid Induced Respiratory Depression
- Pain management
- Oxygen management

"Use of the technology helped drive interventions including airway management and medication optimization for appropriate treatment and avoidance of respiratory depression."

Identifying Patient Deterioration Early Using Contact-Free Continuous Monitoring on an Inpatient Medical-Surgical Unit

Authors: Jan Linderberry MSM, RN and Shelley Derr BSN, RN

DESCRIPTION

Up to 17% of inpatient admissions experience clinical deterioration.¹ Warning signs can often be identified 6 to 8 hours before deterioration events occur.² Early warning scores are used to identify high risk patients, but are not intended to identify deterioration in real time.³ Respiratory and heart rate are the most accurate vital sign deterioration predictors, and accuracy is improved with trending of those vital signs.⁴ Studies have shown contact-free continuous monitoring of respiratory and heart rate on a medical-surgical unit decreases length of stay, ICU length of stay, and code blue rate.⁵

AIM

Implement contact-free continuous monitoring to identify patient deterioration early on a medical-surgical unit.

ACTIONS TAKEN

Contact-free continuous monitoring technology was implemented on a 26-bed inpatient medical-surgical unit at Annot Ogden Medical Center. In addition to providing alerts at customizable upper and lower limits of both respiratory and heart rate, clinicians were able to view trend data for both vital signs.



REFERENCES

1. Ben-art J, Zmilichman E, Adl N, Sorkine P. Contactless respiratory and heart rate monitoring: validation of an innovative tool. J Med Eng Technol. 2010;34(7):8:393-8.
2. Brown H, Terrence J, Vasquez P, Bates DW, Zmilichman E. Continuous monitoring in an inpatient medical-surgical unit: a controlled clinical trial. Am J Med. 2014;127(3):226-32.
3. Prythorch DR, Smith GB, Schmidt PE, Featherstone PL. VIEWS—towards a national early warning score for detecting adult inpatient deterioration. Resuscitation. 2010;81(8):932-7.
4. Churpek MM, Adhikari R, Edelson DP. The value of vital sign trends for detecting clinical deterioration on the wards. Resuscitation. 2016;102:1-5.

SUMMARY OF RESULTS

In this case series, we present seven cases in which contact-free continuous monitoring was used to identify patient deterioration and intervene earlier. Use of the technology helped drive interventions including airway management and medication optimization for appropriate treatment and avoidance of respiratory depression.

Acknowledgments: This study was supported by Heart Rate & Respiratory Rate Monitoring on Critical Care Units - Bed #18 Room, Inc.


CASE TYPE	DEMOGRAPHICS	ADMITTING DIAGNOSIS	RELEVANT HISTORY	CFCM FINDING/ALERT	INTERVENTION/OUTCOME
Airway Management	51-year-old female	Cervical disc, herniation, stenosis	• s/p anterior cervical fusion • esophageal swelling	RR 42	• Contact free continuous monitoring on the bed indicated a high respiratory rate of 42 which triggered the alarm • Brought staff to the bedside and found patient choking but not completely obstructed • Bed showed a one-time burst of elevated respiratory rate due to the patient choking • Ensured patient's airway was patent • Gave patient emotional support – patient expressed relief when staff came in • Prevented patient from hurting herself (she was starting to try to get out of bed) and ensured patient airway
Alcohol Withdrawal Management	22-year-old male	Ethanol withdrawal, Delirium tremens	• Receiving lorazepam for DTs	HR trending up	• Increased lorazepam dose • HR stabilized
Opioid Induced Respiratory Depression	79-year-old male	Total right knee arthroplasty	• Patient recovering from surgery/ anesthesia • Asthma • Periods of apnea	RR 7	• Patient evaluated for sleep apnea • Patient evaluated for naloxone treatment
Opioid Induced Respiratory Depression	69-year-old female	Upper abdominal pain s/p colectomy	• Colectomy for ileocecal adenocarcinoma 1.5 months ago • Guillain-Barre syndrome with residual bilateral leg weakness • Chronic back pain	RR 8	• Patient evaluated for naloxone treatment • Given oxygen • Physician notified
End of Life Care – Pain Management	88-year-old female	Acute pancreatitis	Non-responsive on comfort care	• RR 35 alert • HR/RR both trending up	• Based on trends and physical exam, patient determined to be in pain • Patient given hydromorphone and sublingual atropine • Vitals stabilized • Based on trends and physical exam, pain level improved



Hillrom™

Early Detection of Patient Deterioration Using a Novel Monitoring System¹¹

Early Detection of Patient Deterioration Using Novel Monitoring System



Overview
Patient safety for chronically ill Spinal Cord Injury patients in the hospital is an ongoing challenge. The SCI Center incorporated a novel technology system to detect early patient condition changes and reduce adverse events.

Objective
Describe the EarlySense Monitoring System (ESMS) and discuss its benefits for early detection of patient deterioration in a Spinal Cord Injured (SCI) Medical/Surgical inpatient population.

Participants/Methods
EarlySense provides continuous monitoring of heart and respiratory rate through non-contact technology. Two Spinal Cord Injury med/surg units with similar patient populations were compared during the twelve month period from February 2013 through February 2014 (n=932). This pre implementation data was collected to create a baseline of clinical indicators. During the twelve month period from February 2014 through February 2015 (n=1,150), post implementation data was collected and evaluated in comparison to baseline.

Outcomes
Effect of Safety Initiatives



67%	50%	40%	83%
Medical Response Team (MRT) Activations decreased	Code Blue Activations decreased	ICU Transfers decreased	Mortality following MRT/Code Blue Activations decreased

Enrolled patients were continuously monitored for heart rate (HR) and respiratory rate (RR).

Conclusions
Clinically significant reduction of MRT /Code Blue activations, ICU Transfers, and Mortality was noted among patients on the SCI unit using EarlySense Monitoring System (ESMS).


Poster Presentation
ASCIP Conference
Paralyzed Veterans of America

Clinically significant reduction of MRT/Code Blue activations, ICU Transfers, and Mortality was noted.



Contact-Free, Continuous Monitoring


- ✔ Identify signs of patient deterioration
- ✔ Initiate an early intervention
- ✔ Promote patient safety



Contact-Free, Continuous Monitoring

HOW DOES IT WORK?

When your heart beats
it creates motion.



When you breathe, your lungs expand and contract, creating motion.

The sensor detects cardiac and respiratory motion through the mattress.



The sensor is able to update the HR/RR values twice per second.

All of this information is compiled in an algorithm:

- Creates a running trend of HR/RR
- Filters out other constant motion like an air surface



User Interface Screens



HOME SCREEN

Head Angle 32°

HR 82

RR 18

<30° Limit 30°

Stand Assist

Home, Alarm, Settings, Scale icons

SLEEP SCREEN

Bed Exit On

HR 120

RR 30

Head Angle 32°

Home, Alarm, Settings, Scale icons

HR/RR MONITORING

HR 82

RR 18

40

8

Trends, Display, Suspend buttons

Home, Alarm, Settings, Scale icons

ALERT SCREEN

HR 130

RR 18

82

5

50

6

Low RR Alarm

Respiratory rate exceeds threshold

Time since alarm: 00:27

Suspend button

Home, Alarm, Settings, Scale icons

7-DAY TREND

Apr 25 Apr 26 Apr 28 Apr 30 May 01

HR 81

RR 6

Back, Log, 10m, 1h, 4h, 8h, 24h, 7d buttons


Heart Rate and Respiratory Rate Alerts

When HR/RR exceeds one of the set thresholds, the bed will alert:

VIA LOCAL ALERTS



SafeView®+ Light



Bed Touchscreen Alerts
(turns yellow)




Audible Alert


AND THROUGH NAVICARE® NURSE CALL*



Dome Lights



Status Board



Mobile Device Alerts



SafeView®+ Indicators

SafeView®+ indicator is **WHITE**

- There is no patient in the bed.
- Patient HR/RR can not be read.



SafeView®+ indicator is **GREEN**

- Bed senses patient weight.
- HR/RR is being monitored in safe range.



SafeView®+ indicator is **FLASHING AMBER**

- HR/RR threshold is passed.
- Light turns solid Amber once alarm is silenced.





To learn more about heart rate and respiratory rate monitoring please visit:
<https://www.hillrom.com/Centrella-cfcm>

**CENTRELLA®
SMART+ BED:
THE BED
IS JUST THE
BEGINNING.**



References

1. ECRI Institute. 2019 Top 10 Patient Safety Concerns.
2. Ben-Ari, et al. Contactless respiratory and heart rate monitoring: validation of an innovative tool. J of Med Eng & Tech, Vol. 34, Nos 7-8, Oct-Nov 2010, 393-398.
3. Izrailtyan I, Qiu J, Overdyk F, Erlson M, Gan T. Risk factors for cardiopulmonary and respiratory arrest in medical and surgical hospital patients on opioid analgesics and sedatives.
4. Opioid Safety & Patient Monitoring. The National Coalition to Promote Continuous Monitoring of Patients on Opioids. Nov 2014.
5. Overdyk F, et al. Association of Opioids and Sedatives with Increased Risk of In-Hospital Cardiopulmonary Arrest from an Administrative Database. <https://www.ncbi.nlm.nih.gov/pmc/articles/PMC4767404/>
6. Churpek MM, Adhikari R, Edelson DP. The value of vital sign trends for detecting clinical deterioration on the wards. Resuscitation. 2016;102:1-5.
7. Brown HV, et al. The American Journal of Medicine. 2014; 127:226-232.
8. Sepsis Alliance. Sepsis.org. Sepsis Fact Sheet.
9. Patient Safety Movement; Patient Stories; <https://patientsafetymovement.org/advocacy/patients-and-families/patient-stories/john-lachance/>.
10. Page D MD, et al. Community-, Healthcare- and Hospital-Acquired Severe Sepsis Hospitalizations in the University HealthSystem Consortium. Crit Care Med. 2015 Sept; 43(9): 1945-1951.
11. Early Detection of Patient Deterioration Using Novel Monitoring System. VA Poster. ASCIP Conference.
12. Jones D, et al. Defining clinical deterioration. 2013.
13. Edelson D MD MS FAHA. Agile MD-Hill-Rom Collaboration. Dec 2018.