



**Hillrom™**

**Welch Allyn®  
Connex® Spot Monitor**

# MAKE SMARTER, MORE INFORMED CLINICAL DECISIONS

## Meet the Welch Allyn® Connex® Spot Monitor

It's easy to use and supports comprehensive bedside monitoring—from routine vital signs acquisition to early warning scoring and interval monitoring. And we can integrate with your EMR, giving your care teams immediate access to critical information—when and where they need it.

## DESIGNED WITH PATIENT OUTCOMES IN MIND.

1

Customized features help clinicians identify and respond to signs of patient deterioration sooner

2

Vital signs and other patient health data points are sent directly to the EMR to help reduce errors and omissions

3

Industry-leading technologies support and produce readings you can rely on



**Our guides and features allow you to customize and standardize your workflows.**

Select from six customizable Early Warning Score protocols (eg. MEWS, PEWS, NEWS2) to help detect patient deterioration

Large color touchscreen and single sign on capability simplifies training and operation

Manually enter weight, height, respiration rate, pain level and document modifiers, like body position and O<sub>2</sub> therapy



Spot-check respiration measurement with the Masimo® RRp technology

Our SureBP® technology captures BP readings on inflation in about 15-seconds

Select Masimo®, Nellcor® or Nonin® SpO<sub>2</sub>

Choose Braun ThermoScan® PRO 6000 ear or SureTemp Plus oral/axillary/rectal thermometry

Works with a wide range of cuffs to screen adult, pediatric and neonatal patients

## VISUALLY ASSESSING RESPIRATION RATE? GO DIGITAL.

Digitally spot-check a patient's respiration rate without changing your workflow. Acquire accurate, connected readings at the fingertip helping your team reduce transcription errors and detect signs of deterioration.



## Secure your network and patient health data

If data security is your top concern, we have you covered. When you equip your Connex Spot Monitor with Welch Allyn EMR connectivity solutions, we align with your network security protocols to help you protect your patients and infrastructure.

- FIPS-approved encryption of transmitted data with device and server authentication
- Digitally signed software images and configuration files protect devices from unauthorized changes
- Single sign-on compatibility, including Imprivata, supports quick, secure logins
- Screen lock feature protects idle devices from unauthorized access
- Patient data is automatically deleted after it's sent to the EMR and cleared from the device



DOD RMF



APPROVED





# Hillrom™

## ACCESSORIES

7000-APM	Accessory Power Management Stand
7000-MS3	Classic Mobile Stand
7000-MWS	Mobile Work Surface Stand
7000-GCX	GCX Wall Channel Mount
6000-NC	VSM 6000 Nurse Call Cable
7000-916HS	2D Barcode Scanner
7000-916HSR	2D Barcode Scanner + High Frequency RFID Reader
BATT22	CSM Lithium-Ion Battery
BATT99	APM Stand Lithium-Ion Battery
7000-35	Blood Pressure Hose with FlexiPort® (10 ft)
7000-33	Neonatal Blood Pressure Hose (10 ft)
Extended Warranty#	Extension of the manufacturing warranty
SmartCare Prevention*	Preventive maintenance services
SmartCare Response*	Comprehensive repair services
SmartCare Complete*	Preventive maintenance and comprehensive repair services

# Only available point of sale

\* Service is provided at the customer site and/or at the Hillrom EMEA service center

\* Only available in selected markets.

The Connex Spot Monitor is compatible with a wide range of cleaning products to help you reduce the risk of cross-contamination.



## ACCESSORY POWER MANAGEMENT STAND

7000-APM

Extend your battery life to up to 17 hours of use, helping minimize downtime due to re-charging. Plus, the integrated work surface and built-in accessory/cleaning storage were designed to enhance your productivity.

Want to learn more? Contact your Hillrom account representative today for a demo or visit [hillrom.com/welchallyn](http://hillrom.com/welchallyn)

[hillrom.com](http://hillrom.com)

Navan Business Park, Dublin Road, Navan, Co Meath, C15 AW22, Ireland

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

The Connex Spot Monitor is intended to be used by clinicians and medically qualified personnel for monitoring vital signs of all patients. Medical device (93/42/EEC): Class IIa Notified Body: DQS Medizinprodukte GmbH (CE0297) Manufacturer: Welch Allyn, Inc., 4341 State Street Road, Skaneateles Falls, NY 13153 USA. EC Rep.: Welch Allyn Limited, Navan Business Park, Dublin Road, Navan, County Meath, Republic of Ireland. This medical device is a regulated health product which, pursuant to such regulation bears a CE mark. Welch Allyn recommends that you carefully read the detailed instructions for safe and proper use included in the documents accompanying the medical devices. The personnel of healthcare establishments are responsible for the proper use and maintenance of these medical devices. Welch Allyn reserves the right to make changes without notice in design, specifications and models. The only warranty Welch Allyn makes is the express written warranty extended on the sale its products.

© 2020 Welch Allyn Inc. ALL RIGHTS RESERVED. APR58101-EN-EMEA-R1 04 MAY 2020 ENG - EMEA