

Hillrom™

**Welch Allyn
Spot Vital Signs 4400 Device**

Simplicity, delivered

You work in a fast-paced environment, so you can't afford to slow down when it comes to capturing the basics. Our next-generation vital signs device is here to help. Experience simple, accurate vitals capture with a solution that fits smoothly into your daily routine. The Spot Vital Signs 4400 device offers an efficient way to capture, access and document patient vitals so you can focus on what really matters—spending more time with your patients.

Simple

An intuitive, all-in-one design makes the Spot Vital Signs device easy to use—for caregivers at any level. Add it to your existing workflow without disruption and begin exams quickly from the touchscreen.

Accurate

Trusted for clinical accuracy and consistent measurements,¹ the Spot Vital Signs device can help you make informed treatment decisions at the point of care. The portable design can also help standardize vitals capture throughout your facility and increase productivity for high-volume environments. Plus, blood pressure averaging is now part of the main workflow, enabling you to take and average multiple readings for improved hypertension detection.

Connected

Save time and help reduce data entry errors by sending patient data to the EMR.




Easy to Navigate,
Easy to Use



All-In-One Vitals Capture



EMR Connectivity
via USB Cable





Blood Pressure Averaging with One Tap



Your New Solution for Detecting Hypertension

You may have heard about blood pressure averaging, but do you know how it can help improve hypertension screening?

Relying on a single blood pressure measurement can be misleading. First-time readings can yield higher results than subsequent blood pressure measurements.² Additionally, it's common for variations to occur between readings captured in a series.³ That's why the American Heart Association recommends automated office blood pressure solutions that can capture and average multiple readings during a single visit.²

Blood pressure averaging helps provide a more accurate view of your patient's heart health, so you can detect hypertension with confidence.

EXPERIENCE BLOOD PRESSURE READINGS YOU CAN BELIEVE IN.



ORDERING INFORMATION*

Part Number	Description
44XT	Spot 4400 with NIBP and SureTemp Thermometer
44WT	Spot 4400 with NIBP, SureTemp Thermometer and Nonin SpO ₂
4400-DST	Spot 4400 Desk Stand
4400-MBS	Spot 4400 Mobile Stand
77794-M4400	Green Series 777 Integrated Wall System for Spot 4400
77794-2M4400	Green Series 777 Integrated Wall System with PanOptic Ophthalmoscope for Spot 4400
BATT22	Connex Spot, 2-Cell LI-Ion Battery
106275	Micro USB Cable
REUSE-11	Cuff, WA, Reusable, Adult
REUSE-12	Cuff, WA, Reusable, Large Adult
4500-35	Fast BP Hose with FlexiPort , 10 ft

For more information, please contact your local distributor or Hillrom sales representative.

hillrom.com

Hillrom reserves the right to make changes without notice in design, specifications and models. The only warranty Hillrom makes is the express written warranty extended on the sale or rental of its products.

*See a full list of parts and accessories online at hillrom.com/spot4400.

¹ 2003 std-Alpert BS, Validation of the Welch Allyn Spot Vital Signs blood pressure device according to the ANSI/AAMI SP10: 2002. Accuracy and cost-efficiency successfully combined. Blood Press Monit. 2007 Oct;12(3):345-7.

² Muntner, Paul et al. "Measurement of Blood Pressure in Humans: A Scientific Statement From the American Heart Association." Hypertension 73, no. 5. (May 1, 2019). Doi:10.1161/hyp.0000000000000087.

³ Robert Smith, MD. Blood pressure averaging methodology: Decreasing the risk of misdiagnosing hypertension.