



HillromTM

Welch Allyn Connex ProBP 3400

Accurate, reliable, connected blood pressure readings in the palm of your hand

The Welch Allyn Connex ProBP 3400 puts consistent, accurate blood pressure readings at everyone's fingertips.

- Small, handheld design makes device easy to carry and capture readings on rounds or in triage settings
- SureBP NIBP technology enables fast, accurate 15-second blood pressure readings, improving patient comfort
- Obtains accurate readings and minimizes error codes, even in the presence of many cases of motion or weak pulse
- USB-wired technology enables integration to select EHR/EMR systems, helping improve documentation accuracy and efficiency
- Compatible with Welch Allyn cuff sizes ranging from small child (size 8) to thigh (size 13)
- Lithium-ion battery technology delivers up to 100 readings on a single charge
- Configurable to match your workflow and specific patient needs. Upgradeable so that your initial investment is protected as your needs change



Configured for your facility and your workflow

				
Welch Allyn Connex ProBP 3400	Handheld	Wall Mount	Mobile Stand	Desk Mount
SureBP	34XFHT-B*	34XFHT-B*	34XFST-B*	4602-61 Desk Mount only requires handheld ProBP (not included)
Standard BP	34XXHT-B*	34XXHT-B*	34XXST-B*	

*Power-cord configurations: Replace "B" in part number with designation from chart below.	
North America	B
Europe	2
UK	4
Australia (orange cord)	6
China	C
South Africa	7
Japan	J+

The perfect blood pressure solution for just about any care setting.



Physician Office

With AAMI/BHS-AA-certified accuracy, **SureBP** blood pressure technology that captures readings in just 15 seconds, and the ability to work with a broad range of cuff sizes¹, the **ProBP 3400** is a reliable alternative to manual sphygmomanometers. The ProBP 3400 is also an affordable approach to digitally capture blood pressure measurements and export readings to electronic health record systems.



Hospital

Small enough to carry in your hand, the ProBP 3400 is a convenient alternative to manual sphygmomanometers or larger vital signs devices for capturing reliable blood pressure readings during rounds or in triage settings. A mobile stand featuring a compact footprint and accessory basket enhances portability with an added measure of security and the ability to store extra cuff sizes on board.



Specialty Care

In situations when you need accurate blood pressure readings quickly and easily, look to ProBP 3400 first. Whether in long-term care, sports medicine, dental health, or other applications that call for reliable, repeatable blood pressure readings, the portable and accurate ProBP 3400 is the right choice for you.

Proven Accuracy and Reliability

Using the same ground-breaking **SureBP** technology found in the Welch Allyn Spot Vital Signs LXi, the **ProBP 3400** can capture blood pressure reading in just 15 seconds — nearly twice as fast as other electronic blood pressure devices.

- Meets AAMI standard for blood pressure accuracy
- Achieved “AA” rating by British Hypertension Society data analysis

“...(SureBP technology) offers great advantages in speed and patient comfort...clinicians can have added confidence in this new technology.”*

- **Dr. Bruce S. Alpert**

*Alpert, Bruce S., Clinical evaluation of the Welch Allyn SureBP® algorithm for automated blood pressure measurement. Blood Pressure Monitoring, Vol. 12 No. 4, August 2007.

Connex ProBP 3400 Digital Blood Pressure Device



Connect to your electronic health record system.

The Connex ProBP 3400 comes connectivity-ready, allowing you to integrate it with the electronic health records system you have now or in the future. With hardwired USB connection, the Connex ProBP 3400 communicates with select EHR/EMR systems in both the physician's office and the hospital, helping reduce documentation errors and improve efficiency.



The correct cuff means the correct reading.

Unlike many digital blood pressure devices, the Connex ProBP 3400 works with a very broad range of cuff sizes¹, minimizing the potential for inaccurate readings on patients large or small. Plus you have all the advantages of unique Welch Allyn **FlexiPort** cuffs that enable quick, snap-fit connections, disposable cuffs to reduce cross-contamination risk and facility-standardization benefits.

The ProBP 3400 is not intended for neonatal use.

About Baxter

Every day, millions of patients, caregivers and healthcare providers rely on Baxter's leading portfolio of diagnostic, critical care, kidney care, nutrition, hospital and surgical products used across patient homes, hospitals, physician offices and other sites of care. For 90 years, we've been operating at the critical intersection where innovations that save and sustain lives meet the healthcare providers who make it happen. With products, digital health solutions and therapies available in more than 100 countries, Baxter's employees worldwide are now building upon the company's rich heritage of medical breakthroughs to advance the next generation of transformative healthcare innovations. To learn more, visit www.baxter.com and follow us on Twitter, LinkedIn and Facebook.

For more information, please contact your local distributor or Hillrom sales representative.

hillrom.com

¹ Welch Allyn Motion Tolerance Clinical Study, August 2016

² Alpert, Bruce S., Clinical evaluation of the Welch Allyn **SureBP** algorithm for automated blood pressure measurement. Blood Pressure Monitoring, Vol. 12 No. 4, August 2007.

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.