



Hillrom™



**Welch Allyn SureTemp Plus
692 Electronic Thermometer**

The technology that sets the standard for speed,
accuracy and patient safety with advanced security features

Hillrom is a part of Baxter

Baxter

SURETEMP PLUS 692

WITH THE SPEED AND ACCURACY YOU REQUIRE, AND THE SAFETY AND SECURITY YOU DEMAND

1 EASY

Intuitive user interface, easy-to-read LCD screen and flashing icons, helps select and identify temperature measurements and modes.

2 SAFE

Designed with interchangeable probes and probe wells (blue for oral/axillary and red for rectal) to address patient safety and help reduce risk of cross-contamination.

3 FAST AND ACCURATE

In addition to improved accuracy, the SureTemp Plus 692 is one of the fastest thermometers on the market, providing approximately 4-6 second oral, 10-13 second pediatric axillary and rectal and 12-15 second adult axillary measurements.

4 VERSATILE

The security features make sure the unit stays where it is supposed to be. A location ID on the display and other innovative features help prevent theft.



ONE OF
THE **FASTER**
THERMOMETERS

ORAL
READINGS



4-6 Seconds
approx.

AXILLARY
READINGS



10-13 Seconds
approx.

RECTAL
READINGS



12-15 Seconds
approx.

10-13 Seconds
approx.

SURETEMP PLUS 692

SETTING THE STANDARD OF CARE



CHOICE OF CONFIGURATIONS AND ACCESSORIES:



Handheld



Wall-mounted



Integrated into vital signs monitors



Interchangeable probes and wells with two different lengths



Soft disposable probe covers



In-house calibration checker



Welch Allyn SureTemp Plus 692 Electronic Thermometer

PRODUCT SPECIFICATIONS

Calibration Accuracy	±0.1 °C (±0.2 °F) (Monitor Mode)	
Temperature Predict Time	Oral: Adult Axillary (18 years and older): Pediatric Axillary (17 years and younger): Rectal:	Approximately 4 - 6 seconds Approximately 12 - 15 seconds Approximately 10 - 13 seconds Approximately 10 - 13 seconds
Unit Dimensions	215 mm x 81 mm x 62 mm (8.46" x 3.18" x 2.43")	
Weight	357 g (12 oz)	
Ambient Temperature Range	10 °C (50.0 °F) to 40.0 °C (104 °F)	
Patient Temperature Range	26.7 °C (80.0 °F) to 43.3 °C (110 °F)	
Power Source	(3) 1.5 Vdc AA batteries	
Battery Life	Approximately 6,000 temperatures	

ORDERING INFORMATION

THERMOMETERS

01692-400	SureTemp Plus 692 Electronic Thermometer, Wall Mount, Security System with ID Location Field, Oral Probe and Oral Probe Well with 1.2 m cord
01692-401	SureTemp Plus 692 Electronic Thermometer, Wall Mount, Security System with ID Location Field, Rectal Probe and Rectal Probe Well with 1.2 m cord
01692-500	SureTemp Plus 692 Electronic Thermometer, Wall Mount, Oral Probe and Oral Probe Well with 2.7 m cord (must be locked to the wall)
01692-501	SureTemp Plus 692 Electronic Thermometer, Wall Mount, Rectal Probe and Rectal Probe Well with 2.7 m cord (must be locked to the wall)

ACCESSORIES

21326-0000	Wall Holder
21330-0000	Holder for Spare Probe and Well
06138-000	Calibration Key Assembly, 690/692
02892-000	Probe and Well Kit, 1,2 m, Rectal
02892-100	Probe and Well Kit, 2,7 m, Rectal
02893-000	Probe and Well Kit, 1,2 m, Oral
02893-100	Probe and Well Kit, 2,7 m, Oral
02891-0000	Probe Well, Blue
02891-1000	Probe Well, Red

WARRANTY

S4-M690	Extended Warranty 690/692, 1 year
S4-M690-2	Extended Warranty 690/692, 2 years
S4-M690-5	Extended Warranty 690/692, 2 years

REPLACEMENT PROBE COVERS

05031-101	Case of 1,000 Probe Covers
-----------	----------------------------

For further information about our products and services, please contact your local Hillrom representative or visit our website.

hillrom.com

The Welch Allyn SureTemp Plus 692 thermometer measures patient's temperature in the oral, axillary, and rectal body sites of adult and pediatric patients. Medical device (93/42/EEC): Class IIa. Notified Body: DQS Medizinprodukte GmbH (CE0297) Manufacturer: Welch Allyn, Inc., 4341 State Street Road, Skaneateles Falls, NY 13153 USA. EC Rep.: Welch Allyn Limited, Navan Business Park, Dublin Road, Navan, Co. Meath, C15 AW22, Ireland. This medical device is a regulated health product which, pursuant to such regulation bears a CE mark. Welch Allyn recommends that you carefully read the detailed instructions for safe and proper use included in the documents accompanying the medical devices. The personnel of healthcare establishments are responsible for the proper use and maintenance of these medical devices.

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.