



Hillrom™



Welch Allyn® KleenSpec®
Disposable Vaginal Specula and
Cordless Illumination System

An advancement in women's health
technology that's nothing short of brilliant

Baxter

Hillrom is a part of Baxter

Welch Allyn KleenSpec® Disposable Vaginal Specula and Cordless Illumination System— now in every size



590 SERIES KLEENSPEC® DISPOSABLE VAGINAL SPECULA ENHANCE PATIENT COMFORT

- The extra small speculum allows you to perform a comfortable exam with virginal/pediatric patients as well as postmenopausal and post-hysterectomy women—no more referrals or uncomfortable exams using metal specula or the wrong size.
- Smooth, molded acrylic designed to prevent pinching or binding of tissue for greater patient comfort or the wrong size.

HELP REDUCE CROSS-CONTAMINATION

- Single-use speculum helps to reduce crosscontamination for improved patient safety while eliminating the time and expense of reprocessing metal specs

IMPROVE STAFF PRODUCTIVITY

- Features a wider handle design for improved maneuverability and better balance during exams
- A strong grip and a yoke that's more flexible than low cost substitutes



Look for the distinctive Welch Allyn blue, green, white and yellow colored handles



"We have been using the Welch Allyn cordless Illumination system with disposable vaginal speculums for all of our patient exams since 2008 and couldn't be happier. This light is brighter than any I have ever seen and shows true tissue color. If you forget to charge the light, it is easy to grab a replacement from another room or use a spare. The light is small, easy and quick to use which truly makes us more efficient. I would recommend this system to any clinicians looking to improve the quality of their care."

Shelley Kolton, MD / Downtown Women OB/Gyn Associates / New York City, NY

KLEENSPEC CORDLESS ILLUMINATION SYSTEM

- LED technology as opposed to halogen, for a whiter, brighter light and true tissue color during an exam
- No bulb replacement necessary for a reduced cost of ownership
- Features a cordless design so you'll never have to worry about cords getting in your way, breaking or needing to be cleaned, speeding the exam and giving you a full range of motion
- Push-button on/off; the illumination system easily fits into every size 590 Series speculum
- Newly designed light pipe provides a homogeneous spot and eliminates glare back so all of the light is projected forward
- Uniform spot—no dark or hot spots—for enhanced visualization of the exam area
- Produces much less heat than traditional lamps and has a continuous on time of 80 minutes
- A lithium ion battery provides a recharge time of five hours



**Hillrom™**

Welch Allyn® KleenSpec® Disposable Vaginal Specula and Cordless Illumination System

KLEENSPEC VAGINAL SPECULA—PREMIUM 590 SERIES

590XS	KleenSpec 590 Series Premium Disposable Vaginal Specula, X-Small (24/box; 4 boxes/case)
59000	KleenSpec 590 Series Premium Disposable Vaginal Specula, Small (24/box; 4 boxes/case)
59001	KleenSpec 590 Series Premium Disposable Vaginal Specula, Medium (24/box; 4 boxes/case)
59004	KleenSpec 590 Series Premium Disposable Vaginal Specula, Large (18/box; 4 boxes/case)

KLEENSPEC VAGINAL SPECULA—PREMIUM 590 SERIES WITH SMOKE TUBE

59005	KleenSpec 590 Series Premium Disposable Vaginal Specula with Smoke Tube, Small (12/box; 4 boxes/case)
59006	KleenSpec 590 Series Premium Disposable Vaginal Specula with Smoke Tube, Medium (12/box; 4 boxes/case)

KLEENSPEC CORDLESS ILLUMINATION SYSTEM

80015	Complete KleenSpec Cordless Illuminator with Charging Station, EU and UK plugs (conforms to IEC 60601-1)
80000	KleenSpec Cordless Illuminator
74015	Charging Station for Cordless Illuminator, EU and UK plugs
59010	Disposable Sheath for Cordless Illumination System (100/box; 5 boxes/case)

KLEENSPEC VAGINAL SPECULA—PREMIUM 580 SERIES W/SHEATH

58000S	KleenSpec 580 Series Premium Disposable Vaginal Specula with Sheath, Small (18/box; 4 boxes/case)
58001S	KleenSpec 580 Series Premium Disposable Vaginal Specula with Sheath, Medium (18/box; 4 boxes/case)
58004S	KleenSpec 580 Series Premium Disposable Vaginal Specula with Sheath, Large (12/box; 4 boxes/case)

For further information about our products and services, please contact your local Hillrom representative or visit our website.

hillrom.com

The KleenSpec Disposable Vaginal Specula are used to dilate the vagina and expose the interior of the vagina and exterior of the cervix during pelvic examinations and other gynecological procedures. The KleenSpec Cordless Illuminator provides illumination during pelvic examinations and other general-purpose procedures. Medical Device (EU 2017/745): Class I. Manufacturer: Welch Allyn, Inc., 4341 State Street Road, Skaneateles Falls, NY 13153 USA. EC Rep.: Welch Allyn Limited, Navan Business Park, Dublin Road, Navan, County Meath, C15 AW22, Republic of Ireland. This medical device is a regulated health product which, pursuant to such regulation bears a CE mark. Welch Allyn recommends that you carefully read the detailed instructions for safe and proper use included in the documents accompanying the medical devices. The personnel of healthcare establishments are responsible for the proper use and maintenance of these medical devices.

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

© 2022 Welch Allyn, Inc. ALL RIGHTS RESERVED. APR296101-EN-EMEA-R1 24-OCT-2022 ENGLISH – EMEA