

# IHE Integration Statement

Vendor	Product Name	Version	Date
Welch Allyn	Connex® VM	v2.0	20-Jan-2012

This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:

Profiles Implemented	Actors Implemented	Options Implemented
Consistent Time	Time Client	Secure NTP.
Device Enterprise Communication (2011)	Device Observation Reporter	
Patient Administration Management	Patient Demographics Consumer	Merge
	Patient Encounter Consumer	Inpatient / Outpatient Encounter Management

Internet address for vendor's IHE information: <a href="http://www.welchallyn.com/documents/Patient%20Monitoring/Vital%20Signs%20Capture/Connex%20Data%20Management%20System/IHE%20Integration%20Statement.pdf">http://www.welchallyn.com/documents/Patient%20Monitoring/Vital%20Signs%20Capture/Connex%20Data%20Management%20System/IHE%20Integration%20Statement.pdf</a>		
<b>Links to Standards Conformance Statements for the Implementation</b>		
HL7	<a href="http://www.welchallyn.com/products/en-us/x-18-mu-96-1179856474703">http://www.welchallyn.com/products/en-us/x-18-mu-96-1179856474703</a> .	
DICOM		
<b>Links to general information on IHE</b>		
In North America: <a href="http://www.ihe.net">www.ihe.net</a>	In Europe: <a href="http://www.ihe-europe.net">www.ihe-europe.net</a>	In Japan: <a href="http://www.ihe-j.org">www.ihe-j.org</a>