

# Hear the Difference

*from the leader in diagnostic care.*

FOR VETERINARIANS

## Welch Allyn Veterinary Stethoscopes

For over a decade, we have been focused on providing veterinarians with exceptional performance and value from every Welch Allyn device. Our stethoscopes are no exception.



### Harvey™ DLX™ Stethoscopes

Trumpet brass construction translates to hearing more of the important cardiac, pulmonary, vascular and abdominal sounds without relying on difficult “variable pressure” techniques.

- Available with double- or triple-head rotatable design
- Chrome-plated brass chestpiece and binaurals
- Rotatable binaurals and interchangeable comfort-sealing ear tips
- Includes pediatric diaphragm for use with small animals
- 28" tubing length



#### 5079-315 Accessory Kit

(Includes one pediatric chestpiece, one flat diaphragm, one corrugated diaphragm, one pair soft eartips, one pair hard eartips)



### Harvey™ Elite® Stethoscope

Superior comfort and a stainless steel head that transmits sounds crisply and clearly.

- Double head (bell and diaphragm)
- Rotatable binaurals and interchangeable comfort-sealing ear tips
- Pediatric diaphragm available
- 28" tubing length



#### 5079-232 Accessory Kit

(Includes one pediatric chestpiece, one flat diaphragm, one corrugated diaphragm, one pair soft eartips, one pair hard eartips)

| Product Description | HARVEY DLX   | HARVEY ELITE   |
|---------------------|--|--|
| Product Type        | Double/Triple Head   | Double Head  |
| Head Diameter       | 1.9" (48.3 mm)   | 1.9" (48.3 mm)   |
| Overall Length      | 28" (711 mm)   | 28" (711 mm)   |
| Tubing Length       | 16.5" (419 mm)   | 16.5" (419 mm)   |
| Latex Free          | Yes  | Yes  |
| Chestpiece Material | Brass  | Stainless Steel  |
| Diaphragm Material  | Polycarbonate  | Polycarbonate  |
| Diaphragm Style     | Flat or Corrugated   | Flat or Corrugated   |
| Type of Tubing      | Dual Lumen   | Dual Lumen   |
| Anti-Chill Ring     | No   | Yes  |
| Warranty            | 10-year  | 10-year  |
| Product Numbers     | <ul style="list-style-type: none"> <li>● 5079-321 Triple, Black</li> <li>● 5079-323 Triple, Navy</li> <li>● 5079-325 Double, Black</li> <li>● 5079-327 Double, Navy</li> <li>● 5079-325P Pediatric, Double, Black</li> <li>● 5079-322 Triple, Burgundy</li> <li>● 5079-324 Triple, Forest Green</li> <li>● 5079-326 Double, Burgundy</li> <li>● 5079-328 Double, Forest Green</li> </ul> | <ul style="list-style-type: none"> <li>● 5079-125 Elite, Black</li> <li>● 5079-270 Elite, Burgundy</li> <li>● 5079-271 Elite, Navy</li> <li>● 5079-284 Elite, Forest Green</li> <li>● 5079-125P Pediatric, Elite, Black</li> </ul> |

To learn more about Welch Allyn Stethoscopes contact your local distributor rep, or call 1.800.535.6663.

4341 State Street Road, Skaneateles Falls, NY 13153 USA (p) 800.535.6663 (f) 315.685.3361

© 2012 Welch Allyn MC9131

[www.welchallyn.com/dvm](http://www.welchallyn.com/dvm)

**WelchAllyn®**

Advancing Frontline Care™